



For research use only

ISO9001

100bp DNA Mass Ladder Marker, Ready-to-use

Product	Volume	Cat. No.	Remarks
100bp DNA Mass Ladder Marker, Ready-to-use	500 µl (35 µg)	EBM-1010	For measurement of DNA quantity in agarose gel

Description

The 100bp DNA mass ladder marker is a mixture of double-stranded DNA fragments for determining both the exact size and the amount of DNA fragments after agarose or polyacrylamide gel electrophoresis. The 100bp DNA mass ladder marker consists of DNA fragments ranging in size from 100 to 1,000 bp in 100 bp increments, and additional 1,500 bp fragment. The 100bp DNA mass ladder marker is supplied in a ready-to-use format. This ladder marker can be stained with ethidium bromide or any other known DNA staining methods.

Storage Buffer

- Marker DNA : 35 µg in 0.5 ml of 10 mM Tris-HCl, pH8.0, 1 mM EDTA, 5% Glycerol, 0.005% Bromophenol Blue, and 0.005% Xylene Cyanol

Recommended Storage Condition

- 20°C for 2 year
- 4°C for 6 months
- Room temperature (20-25°C) for 2 months


Usage Information

- Concentration : 700 ng/10 µl (35 µg/500 µl)
- Recommended loading : 5-10 µl (50-100 lanes, ready-to-use)
- Range : 100 – 1,500 bp
- Number of bands : 11

Cautions


- Always use the fresh tip to take out marker solution.
(If you do not, trace amount of contaminated DNases from buffer tank may degrade marker DNA rapidly)
- Don't boil the product.
- Use appropriate % of gels for separation of 100 to 1,500 bp sizes
(1 to 3% agarose gel is recommended)
- Confirm that the concentration of DNA staining dye is optimal before use.
(Breakage or suboptimal concentration of ethidium bromide in gel is a main cause of low estimation of marker concentration or your DNA. 5 ng of DNA should be seen in normal condition)
- Loading volume and concentration should be optimized by gel size, well size, and running length.
- Low sized DNA bands can be gradually disappeared as running is progressing.
(This is because some DNA is getting out from gel to buffer during horizontal electrophoresis, not because the DNA concentration is incorrect. This will be the same for your DNA)

100bp DNA Mass Ladder Marker



10 µl/700 ng/lane ;
1.5% agarose in 0.5x TAE, stained with ethidium bromide

Migration Patterns in Different % of Agarose Gels



Recommended Gel Percentages for Separation of Linear DNA

Agarose Gel, %	Range of Separation, bp	Polyacrylamide Gel, %	Range of Separation, bp
0.5	1,000 - 30,000	3.5	100 - 1,000
0.7	800 - 12,000	5	80 - 500
1	500 - 10,000	8	60 - 400
1.2	400 - 7,000	12	40 - 200
1.4	200 - 4,000	20	5 - 100
2	50 - 2,000		

DNA Size Migration with Sample Loading Dyes

Agarose Concentration, %	Xylene cyanol FF	Bromophenol blue	Orange G
0.7 - 1.7	~4000 bp	~300 bp	~50 bp
2.5 - 3.0	~800 bp	~100 bp	~30 bp

Composition of Gel Electrophoresis Buffers

Buffer	Working Concentration	Stock Concentration (per Liter)
Tris-acetate (TAE)	1x 20 mM Tris-acetate 1 mM EDTA	20x Tris base 96.9 g Glacial acetic acid 22.84 ml 0.5 M EDTA (pH8.0) 40 ml
Tris-borate (TBE)	0.5x 45 mM Tris-borate 1 mM EDTA	10x Tris base 108 g Boric acid 55 g 0.5 M EDTA (pH8.0) 40 ml