



For research use only

ISO9001

## 1Kbp DNA Mass Marker, Ready-to-use

Product	Volume	Cat. No.	Remarks
1Kbp DNA Mass Marker, Ready-to-use	500 µl (41 µg)	EBM-1009	For measurement of DNA quantity in agarose gel

### Description

The 1Kbp DNA mass marker is a mixture of 6 double-stranded DNA fragments ranging in size from 200 to 3,000 bp for determining both the concentration and the size of DNA fragments after agarose gel electrophoresis. The 1Kbp DNA mass marker is supplied in a ready-to-use format. This ladder marker can be stained with ethidium bromide or any other known DNA staining methods.

### Storage Buffer

- Marker DNA : 41 µg in 0.5 ml of 10 mM Tris-HCl, pH8.0, 1 mM EDTA, 5% Glycerol, 0.005% Bromophenol Blue, and 0.005% Xylene Cyanol

### Recommended Storage Condition

- 20 °C for 2 year
- 4 °C for 6 months
- Room temperature (20-25 °C) for 2 months


### Usage Information

- Concentration : 820 ng/10 µl (41 µg/500 µl)
- Recommended loading : 5-10 µl (50-100 lanes, ready-to-use)
- Range : 200 – 3,000 bp
- Number of bands : 6

### Cautions


- Always use the fresh tip to take out marker solution.  
(If you do not, trace amount of contaminated DNases from buffer tank may degrade marker DNA rapidly)
- Don't boil the product.
- Use appropriate % of gels for separation of 200 to 3,000 bp sizes  
(0.7 to 1.5% agarose gel is recommended)
- Confirm that the concentration of DNA staining dye is optimal before use.  
(Breakage or suboptimal concentration of ethidium bromide in gel is a main cause of low estimation of marker concentration or your DNA. 5 ng of DNA should be seen in normal condition)
- Loading volume and concentration should be optimized by gel size, well size, and running length.
- Low sized DNA bands can be gradually disappeared as running is progressing.  
(This is because some DNA is getting out from gel to buffer during horizontal electrophoresis, not because the DNA concentration is incorrect. This will be the same for your DNA)

### 1Kbp DNA Mass Ladder Marker



10 µl/820 ng/lane ;  
1.5% agarose in 0.5x TAE, stained with ethidium bromide

### Migration Patterns in Different % of Agarose Gels



### Recommended Gel Percentages for Separation of Linear DNA

Agarose Gel, %	Range of Separation, bp	Polyacrylamide Gel, %	Range of Separation, bp
0.5	1,000 - 30,000	3.5	100 - 1,000
0.7	800 - 12,000	5	80 - 500
1	500 - 10,000	8	60 - 400
1.2	400 - 7,000	12	40 - 200
1.4	200 - 4,000	20	5 - 100
2	50 - 2,000		

### DNA Size Migration with Sample Loading Dyes

Agarose Concentration, %	Xylene cyanol FF	Bromophenol blue	Orange G
0.7 - 1.7	~4000 bp	~300 bp	~50 bp
2.5 - 3.0	~800 bp	~100 bp	~30 bp

### Composition of Gel Electrophoresis Buffers

Buffer	Working Concentration	Stock Concentration (per Liter)
Tris-acetate (TAE)	1x 20 mM Tris-acetate 1 mM EDTA	20x Tris base 96.9 g Glacial acetic acid 22.84 ml 0.5 M EDTA (pH8.0) 40 ml
Tris-borate (TBE)	0.5x 45 mM Tris-borate 1 mM EDTA	10x Tris base 108 g Boric acid 55 g 0.5 M EDTA (pH8.0) 40 ml