



For research use only

ISO9001

rTaq Plus 5x PCR PreMix

Product Description	Quantity	Cat. No.	Remarks
rTaq Plus 5x PCR PreMix	96 reactions	EBT-1321	4 µl/20 µl reaction
	5 x 96 reactions		

Description

rTaq Plus 5x PCR PreMix is a ready-to-use premix pre-aliquoted in PCR tubes containing all the components essential for a PCR and agarose gel electrophoresis (DNA polymerase, dNTP, reaction buffer, glycerol, bromophenol blue, and stabilizer). PCR can be performed simply by adding primer pair, template and PCR grade distilled water.

As PCR PreMix is supplied as a 5x concentration format in 4 µl volume within PCR tubes, users should adjust a final reaction volume to 20 µl by adding primers, template DNAs and distilled water. One unit of rTaq Plus DNA Polymerase is contained in 4 µl of 5x PCR PreMix.

rTaq Plus DNA Polymerase is a thermostable DNA polymerase that designed for a reliable amplification of long and complex targets with a robust yield in an extremely short extension time. Due to its high elongation rate (about 100 bp/sec) and high processivity (about 250 bases), long sized products can be generated within relatively shorter running time compared to any other commercially available enzymes (5 kbp within 60 sec extension).

Recombinant rTaq Plus DNA Polymerase is purified from a modified recombinant *E.coli* strain. The enzyme catalyzes the incorporation of nucleotides into duplex DNA in the 5'→3' direction in the presence of Mg²⁺ at 70–80°C. rTaq Plus DNA Polymerase has no detectable exonuclease activity (5'→3' or 3'→5') and generates a 3'-dA overhang, suitable for a TA cloning of PCR products.

rTaq Plus DNA Polymerase is free from any nucleic acid which can be contaminated from expression host during purification process. That provides robust conditions for amplification of gram-negative bacterial genes, and profiling of these.

QC tests

Performance tests, genomic DNA contamination test, confirmation test for the absence of endo and exonucleases.

Storage Condition

Store at -20°C

rTaq Plus 5x PCR PreMix is stable for at least 2 years at recommended storage condition.

Standard Protocol

1. Prepare 20 µl PCR solution as follows:

5x PreMix in PCR tube :	4 µl
Primer (10 pmol/µl) :	0.5 µl each
Template :	0.1–10 ng
PCR grade distilled water :	– µl

Adjust final vol. to 20 µl with PCR grade distilled water

2. Set PCR cycling as follows :

Initial denature at 95°C : 3 min

		1–2 kbp	3–4 kbp	5–10 kbp
Denature	95°C	10 sec	10 sec	20 sec
Anneal	Tm–4°C	10 sec	10 sec	20 sec
Extend	72°C	10 sec	20 sec	10 sec/kbp

* Set 25–35 PCR cycles for effective amplification

* You can also use two step cycle amplification (denaturation at 95°C and annealing/extension at 68°C)

3. You can analyze PCR products by direct loading into agarose gel because PCR PreMix contains glycerol and bromophenol blue (blue color) essential for a gel loading.

Trouble-Shooting

1. No products

– Confirm your template is intact : Try another reaction with a result assured primer pair and templates

2. Smear bands or smeared background

– Reduce template concentration : High concentration of template can lead to smearing of PCR products. Generally, 1–10 ng of plasmid DNA and 10–100 ng of genomic DNA are working well.
– Increase annealing temperature
– Set up a reaction mix on ice

3. Non-specific bands

– Increase annealing temperature
– Consider using PCR additives, like 1–2% DMSO or 0.5–1 M betaine
– Confirm specificity of your primers

4. Low yield

– Increase PCR cycle number
– Be sure appropriate concentration of your template is added