



For research use only

ISO9001

Bst DNA Polymerase, exo-

| 제 품 명 | 용 량 | Cat. No. | 비 고 |
|--------------------------|------------|----------|-------------------|
| Bst DNA Polymerase, exo- | 6,000 unit | EBT-1701 | 120 unit/ μ l |

Description

Bst DNA Polymerase exo- is a recombinant DNA polymerase from *Bacillus stearothermophilus*. Bst DNA Polymerase is purified from an *E.coli* strain carrying a plasmid with the cloned gene encoding large fragment of DNA polymerase. It contains the 5'→3' polymerase activity, but lacks 5'→3' exonuclease activity. For this reason, Bst DNA Polymerase can be effectively used in chain displacement reaction at high temperature, 65°C.

Applications

- Isothermal amplification (LAMP)
- Applications requiring strand-displacement DNA synthesis
- DNA sequencing through high GC regions
- Rapid sequencing from nanogram amounts of DNA template

Storage Buffer

120 unit/ μ l in 10 mM Tris-HCl, pH8.0, 50 mM KCl, 0.1 mM EDTA, 1 mM DTT, 0.1% Triton X-100, and 50% Glycerol.

10x Reaction Buffer w/ MgSO₄

200 mM Tris-HCl, pH8.8, 100 mM (NH₄)₂SO₄, 700 mM KCl, 30 mM MgSO₄, 0.1% Tween20.

Unit Definition

One unit of enzyme incorporates 10 nmole of dNTP into acid insoluble material in 30 min at 65°C.

QC Tests

Activity, SDS-PAGE purity, performance tests, genomic DNA contamination test, confirmation test for the absence of endo and exonucleases.

Storage Condition

Store at -20°C.

Standard Protocol

- Reaction condition
ex. (LAMP / Loop-mediated isothermal amplification)

| | |
|------------------------------|-------------------------|
| PCR grade distilled water | - μ l |
| 10x Bst reaction buffer | 2 μ l |
| 10 mM dNTP mix (2.5 mM each) | 2 μ l |
| FIP/BIP Primers (20X) | 1 μ l (1.6 μ M) |
| F3/B3 Primers (20X) | 1 μ l (0.2 μ M) |
| LF/LB Primers (20X) | 1 μ l (0.4 μ M) |
| Template DNA | - μ l |
| Bst DNA Polymerase, exo- | 0.2 - 1 μ l |

Adjust final vol. to 20 μ l with PCR grade distilled water

- Incubate reaction at 65°C for 60 min.